

# EMSCUPLT Instruction



When operating this product, please read the manual of this instrument carefully.

## Contents

1. Introduction to the instrument.....	1
1.1 Principle of the instrument ... ..	1
1.2 Instrument parameters.....	1
1.3 Appearance of the instrument.....	1
2. Precautions.....	3
2.1 Precautions before use.....	3
2.2 Contraindications.....	4
2.3 Precautions for customers during and after use.....	4
3. Operation process of the instrument interface.....	5
3.1 Start interface operation flow.....	5
3.2 Handle separate working interface operation.....	8
3.3 Setting interface.....	9
4. Maintenance and repair.....	13
4.1 Maintenance ... ..	13
4.2 Maintenance.....	13
5. Common faults and troubleshooting methods .....	14
6. Responsibilities of the manufacturer .....	15
7. Storage and transportation.....	15
8. Term of use .....	16
9. Parts List .....	17



### 1.1 Principle of the instrument:


Use HI-EMT (High Energy Focused Electromagnetic Wave) technology to make muscles expand and contract continuously, and perform extreme training to deepen the internal structure of muscles, that is, the growth of myofibrils (muscle enlargement) and the production of new protein chains and muscle fibers ( Muscle hyperplasia) to train and increase muscle density and volume. The 100% limit muscle contraction of HI-EMT technology can trigger a large amount of fat decomposition. Fatty acids are broken down from triglycerides and accumulated in large amounts in fat cells. Fatty acids are too high, causing fat cells to apoptosis, which will be excreted from the body through normal metabolism within a few weeks.

### 1.2 Instrument parameters

Name	Parameters
Input power	220V 50HZ
Output power	2800W
Frequency	5-100Hz
Screen Size	17.3 inches
Machine Color	Grey
Handpiece Life	9000 minutes

### 1.3 Appearance of the instrument: front view of the instrument





Note: After the machine is installed, please insert the equipped IC and turn it on for operation.

## 2. Precautions

### 2.1 Precautions before use

- Do not have a full stomach during use, and at least one hour after a meal.
- Before using the instrument, please remove any metal jewelry worn on your body.

- Before starting the operation, first touch the probe vertically to your skin, and then turn on the main unit to change the operation mode  
Adjust the intensity to the lowest level, and then increase the intensity sequentially according to the user's skin bearing capacity;
- Do not eat or drink during the operation, and try to eat one hour later after the end.

## 2.2 Contraindications


- People with wounds on the skin, or avoid the wounds;
- Heart disease patients (such as those with arrhythmia), especially those with pacemakers;
- Cancer patients;
- Infectious disease patients;
- Pregnant women, breastfeeding, menstrual period;
- Major surgery patients, diabetics;
- People with any metal in the body;
- People with cerebrovascular disease (stroke patients)

## 2.3 Precautions for customers during and after use


- After use, the user should replenish water in time;
- Do not overeating, arrange eating and sleeping reasonably;
- Be sure to strengthen moisturizing and sun protection after use;
- Avoid spicy and irritating food;
- Please refrain from hot springs and saunas for three days;

### 3. Operation process of the instrument interface

#### 3.1 Start interface




Note: 1. The operation interface of Mode 1 is as follows (before operation, please take out the handle from the handle holder, and fix the handle on the part required by the customer with a belt)




2. Enter the mode 1 interface to adjust the handle operation time

according to customer needs, and then click the start/stop button to proceed to the next step.




3. Choose according to different genders, ages, and parts (are recommended to choose abdomen) after selecting, click OK to enter the next operation interface.






4. Adjust the corresponding intensity according to different customer needs to start the operation. If the customer requests to pause, you can click the start/pause button to pause the work. After the work is completed, you can directly click the return button to return.




Note: 1. The operation interface of mode 2 is as follows (before operation, please take out the handle from the handle holder, and use a belt to fix the handle on the part required by the customer)




4. Adjust the corresponding intensity/phase frequency/time according to different customer needs. Click the start/pause button to enter the operation work. If the customer requests to pause, you can click the start/pause button to pause the work. After the work is completed, you can directly click the return button return.

### 3.2 Handle separate working interface operation



1. The A/B handle can work alone. If you need to stop the B handle when the handle is working, you can click the mark of the B handle on the screen (as shown in the figure above, the B handle has stopped working, and the A handle is working again) if you need the B handle Restart work, just light up the B handle logo.






2. The above picture shows that the B handle is working, and the A handle stops working.



Note: 1. Align the handle with the operating part and tie it, and then press the start button on the screen to start the operation;



2. After use, gently place the handle on the handle holder;

### 3.3 setting interface









Error code: 1 The handle is not inserted, 2 The service life is exhausted,  
3 The output voltage is too low, 4 The output voltage is too high,  
5 The input AC voltage is too low, 6. The input AC voltage is too high,  
7. The IC card is not inserted 8. The IC card module Communication  
error 9. The IC card is wrong with the handle

#### 4. Maintenance and repair

##### 4.1 Maintenance

4.1.1 Daily inspections Managers must regularly inspect the instruments and components. Whether the power socket is reliably connected to the ground. Please check whether the power cord is damaged before using the instrument, and stop using the instrument if it is damaged.

4.1.2 Cleaning The instrument should be cleaned regularly to ensure the sanitation and good operation of the equipment. You can wipe the housing of the device with a damp cloth and use soapy water or a neutral detergent when the instrument is powered off.

4.1.3 Replacement of the fuse When replacing the fuse, in order to prevent electric shock, please turn off the power and unplug the power plug.

4.2 After repairing the instrument, you need to fill in the warranty card and send it to the manufacturer or agent to have the warranty effect; from the date of issuance of the purchase of the instrument, the product quality problem is guaranteed by the company for one year. According to

the warranty card warranty; the warranty scope here is limited to the instrument host, excluding accessories and other consumables that need to be replaced regularly). After the warranty period has expired, please contact the manufacturer directly to provide paid after-sales service. To repair this machine, please contact our company's business department or dealer, and do not repair it without authorization. If you open the package without permission, our company will not provide warranty.

#### 5. Common faults and troubleshooting methods

Failure phenomenon	Causes	Exclusions and solutions
Device not working	<ul style="list-style-type: none"> <li>● The power plug is not inserted into the socket</li> <li>● The power switch is off</li> <li>● Power fuse failure</li> </ul>	<ul style="list-style-type: none"> <li>● Check power</li> <li>● Check the power switch</li> <li>● Replace the fuse</li> </ul>
The handle does not work	<ul style="list-style-type: none"> <li>● Is the number of probes used up?</li> <li>● The working head is damaged</li> <li>● The interface on the back of the</li> </ul>	<ul style="list-style-type: none"> <li>● Check the number of probes</li> <li>● Replace working head</li> <li>● re-install</li> </ul>



	instrument is not installed correctly	
--	--	--

## 6. The responsibility of the manufacturer

The products produced should meet the product standards. If there is any non-compliance, the company is responsible for replacement and return.

The company shall bear the corresponding responsibility for personal and property damage caused by product manufacturing defects, but the company shall not be liable for the following situations: The purchaser and operating user do not use the instructions in accordance with the operating specifications, and do not follow the safety precautions and warnings. Instructions used.

Violation of product transportation, installation, use, maintenance, storage, etc., causing product, personal, or property damage.

The company implements life-long maintenance and technical services for the products. Any manufacturing quality problems that occur within the specified quality guarantee period will be repaired for users free of charge.

When accessories need to be replaced, original accessories must be purchased from our company, otherwise the adverse consequences will be borne by the user.

## 7. Storage and transportation management

Storage: The instrument should be stored in the original packaging box, placed in a well-ventilated room, the packaging box should be padded, the ambient temperature is  $-10^{\circ}\text{C} \sim 40^{\circ}\text{C}$ , the relative humidity is not more than 93%, and it is not allowed to be placed in harmful gases.

Surrounding flammable, explosive materials and corrosive gases.

Transportation: The instrument in the packaged condition is suitable for road, railway, aviation and waterway transportation. During loading and unloading and transportation, it should be protected from severe vibration and shock, and should not be exposed to moisture, and cannot be mixed with flammable and corrosive substances. , The specific requirements are in accordance with the order contract.

## 8. Use period

The warranty period of this instrument is 1 year (Note: Man-made damage is excluded).

The warranty period for accessories is 3 months (Note: Man-made damage is excluded).

If the instrument is beyond the useful life, it needs to be tested and approved by the engineer, and the performance of the instrument needs to be tested by the engineer regularly in the subsequent use.

Discarded instruments and accessories that have exceeded their service life should be disposed of in accordance with local regulations.

## 9.Parts list

No.	Name	QTY	Picture
1	Main Machine	1pc	
2	Handle	2pcs	
3	Bandage 8*80CM 8*90CM 8*100CM	3pcs	
4	Power wire	1pc	
5	M3/M4 wrench	2pc	
6	Instruction	1pc	<p>EMSCUPLT Instruction</p>  <p><small>When opening the product, please read the manual of this instrument carefully.</small></p>